Beta Bionics[®]

iLet[®] Flex
Infusion Set

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en

Do not use if package is damaged and consult instructions for use

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R_x
ONLY

en

Prescription required, does not apply in all geographies

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Distribuidor

INSTRUCTIONS FOR USE

Read all instructions carefully before using the iLet Flex infusion set.

INDICATIONS FOR USE:

The iLet Flex infusion set is indicated for subcutaneous infusion of medication administered by an external pump. The infusion set is indicated for single use.

DESCRIPTION

The iLet Flex infusion set has a 90-degree soft cannula. It is delivered ready to use in a pre-loaded insertion device with automatic needle retraction. It is provided sterile and non-pyrogenic.

INTENDED USE

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes mellitus. The infusion set has been tested to be compatible with insulin approved for subcutaneous infusion.

CONTRAINDICATIONS

This infusion set is indicated for subcutaneous use only. Do not use the infusion set for intravenous infusion. Do not use the infusion set with blood or blood products.

GENERAL INSTRUCTIONS

- The infusion set may be worn for a maximum of three days, or per a healthcare professional's instructions.
- If the infusion set is being used for the first time,

do the first set-up in the presence of a healthcare professional.

WARNINGS

- Do not use the infusion set if the package is opened or damaged. An opened or damaged package may be contaminated. This can cause an infection.
- Do not use the infusion set if the disconnect cover is removed. Use a new infusion set instead. One purpose of the disconnect cover is as a safety feature against accidental firing of the insertion device which could cause injury.
- Do not change the infusion set just before bedtime unless blood glucose can be checked one to three hours after insertion.
- Do not put alcohol, disinfectants, perfumes, deodorants, cosmetics, or other substances with solvents on the infusion set. These substances may damage the integrity of the infusion set.
- Never point the insertion device toward any body part where insertion is not desired.
- Always rotate the insertion site when changing the infusion set. Reusing an insertion site too often may cause scarring and unpredictable insulin delivery. Review the pump user guide on how to rotate insertion sites.
- Check the insertion site often through the clear window. Improper insertion and maintenance of the insertion site can cause

inaccurate insulin delivery, infection, or site irritation. Replace the infusion set at a new site if the soft cannula is not properly inserted or the adhesive tape becomes loose or displaced.

- Check the cannula housing and tubing for blood. Blood may cause insufficient insulin delivery. This may result in high blood glucose. Replace the infusion set at a new site if there is blood.
- Fill the tubing completely with insulin before insertion. Do not leave air in the tubing. An uncontrolled amount of insulin may be delivered if air is present.
- Check the tubing for clogs or leaks if a blood glucose reading is high. Clogs or leaks may restrict insulin delivery and result in high blood glucose. Replace the infusion set even if a clog or leak is suspected, but not found.
- Do not attempt to clear air or an occlusion in the tubing while it is connected to the body. An uncontrolled amount of insulin may be delivered. This can cause high or low blood glucose. Disconnect the tubing before adjusting it.
- Do not reuse the infusion set. Reuse of the infusion set may cause damage to the cannula or needle and may lead to infection, site irritation, or unpredictable insulin delivery.
- For no delivery alarms, refer to the pump user guide.

PRECAUTIONS

- Do not use another type of infusion set without consulting a healthcare professional for correct handling. A healthcare professional should always be consulted when choosing an infusion set.
- Choose insertion sites as recommended by a healthcare professional. The choice of site depends on treatment and patient-specific factors such as body composition and physical activity level.
- Check that the insertion site is free of skin irritation, such as redness, scar tissue, or bleeding. Do not insert the infusion set into muscle or over bone. Doing so can cause pain or may damage the infusion set. If the infusion set is improperly inserted, replace the infusion set and insert at a new site.
- If needed, remove body hair around the insertion site to ensure that the adhesive sticks to the skin.
- Do not use an insertion site located under a belt or waistband, or where the site is constrained by clothing or accessories. Insertion at these locations can cause disconnection of the infusion set and interruption of insulin delivery, and lead to hyperglycemia.
- Do not reposition the infusion set on the body after the adhesive is placed on the skin. Repositioning the infusion set may damage the adhesive. Replace the infusion set if the adhesive is damaged.
- Always check blood glucose one to three hours after inserting a new infusion set. This is to

confirm accurate insulin delivery. High or low blood glucose can result if insulin delivery is not accurate.

- Check blood glucose several times throughout the day or as recommended by a healthcare professional.
- Consult a healthcare professional for how to correct for missed insulin and for how much time the pump can remain disconnected.

STORAGE AND DISPOSAL

- Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight, or in high humidity.
- Store and handle insulin according to the manufacturer's instructions.
- Dispose of the insertion device in a proper sharps container according to local laws.
- Dispose of a used infusion set according to local regulations for biohazardous waste.

SERIOUS INCIDENT

- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

ADDITIONAL INFORMATION FOR THE HEALTHCARE PROFESSIONAL

- Consider treatment and patient-specific factors when choosing an insertion site.
- The length of the soft cannula should be based on treatment and patient factors such as physiology and activity level. If the soft cannula is too long, insertion pain or insertion into the bone


or muscle may occur. If the soft cannula is too short, leakage or irritation at the insertion site may occur. In both cases, the infusion must be stopped and a new soft cannula with a different length must be inserted at a new insertion site.

- Inform the patient that the selection of an insertion site is dependent on cannula length.


INSTRUCTIONS FOR USE

- Read all instructions carefully before using the iLet Flex infusion set.
- Consult the pump user guide for important information about pump therapy prior to connecting the infusion set. This information includes connection and filling procedures, possible errors, and potential risks related to pump therapy.
- Follow good hygienic procedures. If this infusion set is being used for the first time, perform the first set-up in the presence of a healthcare professional.


THE PRODUCT (1a)

INSERTION DEVICE & CANNULA HOUSING (1b)

TUBING (1c)

1 COMPONENTS

[A] Insertion device

[B] Safeguard (Optional: transparent cover for cannula housing)

[C] Activation button (for insertion of the infusion set)

[D] Indication for tubing direction

[E] Adhesive tape

[F] Paper

[G] Insertion site window

[H] Cannula housing

[I] Soft cannula

[J] Tubing

[K] Site Connector

[L] Pump connector (Luer lock)

[M] White cover for site connector

INSERTING

Unfold the front cover of this booklet to view images that are keyed to the steps below.

- 2 Wash hands with soap and water.
- 3 Select a recommended insertion site (shown in gray) as indicated by a healthcare professional.
- 4 Clean the insertion site with a disinfectant as directed by a healthcare professional. Allow to air dry before inserting the infusion set. Remove hair around the insertion site to ensure the adhesive sticks to the skin if needed.
- 5 Open the package and remove the paper from the tubing.
- 6 Follow the steps to prepare your pump as instructed in the pump user guide. Twist the pump connector (Luer lock) and connector together until tight.
- 7 Remove the white cap from the site connector by (1) squeezing the sides of the site connector then (2) pulling the white cap off. Keep the white cap for later use.
- 8 Connect the pump connector (Luer lock) to the pump as instructed in the pump user guide. Completely fill the tubing with insulin as instructed by the pump user guide. Do not leave any air in the tubing. The tubing is completely filled when drops of insulin are seen coming out of the site connector.
- 9 Remove the paper backing from the adhesive. Be careful not to touch the adhesive.
- 10 Remove the disconnect cover from the insertion device by (1) gently squeezing the sides of the disconnect cover then (2) pulling

it away from the insertion device. Keep the disconnect cover for later use.

- 11** The raised arrows indicate the tubing direction. Stretch the skin until smooth. Then press the insertion device against the skin. Press the top button down completely to insert the infusion set. Gently and carefully remove the insertion device.
- 12** Gently massage the adhesive onto the skin with one finger. Replace the infusion set if the adhesive does not stick to the skin.
- 13** Gently hold the cannula housing steady with one finger. Then push the site connector straight into the cannula housing until it clicks.

Fill the soft cannula with insulin:

6 mm: 0.6 units (0.006 ml)

9 mm: 0.6 units (0.006 ml)

Dispose of the insertion device in an appropriate sharps container according to local laws.

Keep the disconnect cover and the white cap for use when the infusion set is disconnected.

DISCONNECTING

- A** Wash hands with soap and water.
- B** Gently hold the cannula housing steady with one finger. Then (1) squeeze the sides of the site connector and (2) pull the site connector out of the cannula housing.
- C** Put the white cap on the site connector.
- D** Put the disconnect cover on the cannula housing. Then push the cover into the cannula housing until it clicks.

RECONNECTING

- E** Wash hands with soap and water.
- F** Use one finger to gently hold the cannula housing steady. Then (1) squeeze the sides of the disconnect cover and (2) pull the disconnect cover out of the cannula housing.
- G** Remove the white cap from the site connector by (1) squeezing the sides of the site connector then (2) pulling the white cap off. Ensure that there is no air in the tubing. If there are air bubbles in the tubing: Fill the tubing with insulin as instructed in the pump user guide. Do not leave any air in the tubing.
- H** Use one finger to gently hold the cannula housing steady. Then, push the site connector into the cannula housing until it clicks.

REMOVING

- I** Wash hands with soap and water.
- J** Use one finger to gently hold the cannula housing steady. Then (1) squeeze the sides the site connector and (2) pull the site connector out of the cannula housing.
- K** For specific instructions to disconnect tubing from the pump, refer to the pump user guide.
- L** Carefully lift the adhesive around the cannula housing. Then, pull the soft cannula out of the skin.

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Beta Bionics
11 Hughes, Irvine, CA
92618, USA
855-745-3800
www.betabionics.com

© 2025 Beta Bionics, Inc. Beta Bionics® and iLet® are registered trademarks of Beta Bionics, Inc. Bionic Pancreas™ and iLet Bionic Pancreas™ are common law trademarks of Beta Bionics, Inc. All rights reserved. All other trademarks are the property of their respective owners. The use of third-party trademarks does not constitute an endorsement or imply a relationship or other affiliation.